

Small Animal Induction Chambers

An excellent addition to conventional inhalation anesthesia procedures, the Small Animal Induction Chambers are available in four sizes to accommodate cats, small dogs, birds, mice and other small patients.

- Available in 45 liter, 18 liter, 5.25 liter and 1.5 liter sizes
- Constructed of .22" (5.6mm) impact resistant acrylic
- Inlet port accepts standard 22mm breathing tube or 15mm connector
- Outlet port accepts standard 19mm scavenging tube or 15mm connector
- Sliding, self-locking top has stops to prevent unintentional removal of lid
- Scratch-resistant feet
- Includes a hook-up kit to connect vaporizer to inlet port

