

Cutanplast® Sterile haemostatic gelatin




Cutanplast® is a reabsorbable sterile gelatin sponge with a marked haemostatic effect: the porous surface of the gelatin induces the rapid rupture of the blood plaques with the consequent activation of the enzymes cascade which leads to natural coagulation.




Cutanplast® Characteristics

- Completely absorbed if left in-situ.
- Non-toxic.
- Does not cause allergic reactions.
- Has a neutral pH.
- Absorbs a quantity of blood corresponding to approximately 50 times its own weight.
- May be applied dry.
- May be soaked with sterile physiological solution or with antibiotics.
- May be easily reduced in size to any required format.
- Reduces the time of the surgical operation.
- Stops the haemorrhage.
- Is reabsorbed within a number of weeks.
- Easy to use.
- Safe.
- Reliable.


Examples



1 *Neurosurgery*




4 *Chest Surgery*




6 *Parenchymal Surgery*



of where Cutanplast[®] is utilised:




2 *Plastic Surgery*




3 *Otolaryngological Surgery*



5 *Vascular Surgery*



7 *Thyroid Surgery*




8 *Anal Surgery*

TECHNICAL DATA



“Standard”

Dimensions of the sponge mm 70x50x10;
every sponge is individually packaged in
double envelope, sterile.
Packaging boxes of 2 or 20 units.




“Special”

Dimensions of the sponge mm 70x50x1;
every sponge is individually packaged in
double envelope, sterile.
Packaging boxes of 20 units.



“Film”

Dimensions of the sponge mm
200x70x0.5;
every sponge is individually packaged in
double envelope, sterile.
Packaging boxes of 20 units.



“Anal”

Dimensions of the sponge mm 80x Ø 30;
every sponge is individually packaged in
a plastic container, sterile.
Packaging boxes of 5 or 20 units.



“Dental”

Dimensions of the sponge mm 10x10x10;
blister packaging, sterile.
Packaging boxes of 24 units.



Mascia Brunelli S.p.A.

V.le Monza, 272 - 20128 Milano

Tel. ++39 - 02 - 25.20.91 - Fax ++39 - 02 - 25.76.428 - Telex 330838 MABIOS I

Internet: www.masciabrunelli.it - e-mail: mktg@masciabrunelli.it

CE 0483