

Broncho-/ Gastroscopy - LED-Videoendoscope Eickview 70

Art.-No.: 306010



Technical Data:

- Length: 70 cm
 - Diameter: 6,0 mm
 - Working Channel: $\varnothing = 2,0$ mm
 - Depth of field: 7-60 mm
 - Direction of view: 0°
 - Field of view: 90°
 - Deflection: 180° up, 160° down
 - Data format: JPG, AVI
 - Integrated LED light in distal end, power supply comes from monitor
-
- Comes with 3,5" monitor for image and video recording
 - Archiving via SD-card
 - Operation time of monitor: 3 h, comes with charger
 - USB-connection for image and video transfer
 - Connectable to separate monitor via cable
 - Scope is fully immersible