



PET-VET BIOMED

## OXYVET III OXYGEN CONCENTRATOR

An oxygen concentrator draws in room air, separates the oxygen from other gases in the air and delivers the oxygen at high concentrations. This device provides up to 5 l per minute of constant flow oxygen.

- ▶ Flowrate: 0.5 – 5 l/min
- ▶ Oxygen concentration: 93 % ± 3 %
- ▶ Very quiet, only 40 dbA (50 Hz)
- ▶ Low power consumption, 290 watts only
- ▶ Dimensions (in cm): L 34 x W 30 x H 62
- ▶ Weight: 16.3 kg

**3 YEAR WARRANTY**

**213175**

Oxygen Concentrator Bubble Humidifier with Bottle

- ▶ Humidifies oxygen prior to patient delivery

**213176**



213175

**PET-V**

32 Shaham St.  
Petach Tikva, Israel  
Tel: 972.3.9277290  
Fax: 972.3.9277278  
info@petvetBiomed.com  
www.petvetBiomed.com

8 Hamelacha St.  
Entrance A  
2nd floor, Suite 206  
South Industrial Area  
Netanya

רח' המלאכה 8  
כניסה א', קומה ב'  
חדר 206  
אזור תעשייה דרומי  
נתניה

רח' שחם 32, קרית מטלון  
ת.ד. 7803, פתח תקוה 4951727  
טל.: 03.9277290  
פקס: 03.9277278  
info@petvetBiomed.com  
www.petvetBiomed.com