



## PET-VET BIOMED

### 194022- Lateral Suture Set – Basic

Cranial cruciate ligament rupture is one of the most common reasons for hind limb lameness in dogs. Extra-capsular suture stabilisation with monofilament nylon is a well-established repair technique for this condition. This surgical stabilisation set provides all the necessary equipment for the veterinary surgeon to perform this procedure.

Includes:

#### Crimping Forceps

#### 2 x CCL Suture System (Swaged Needle)

- 30 mm fabella needle
- 50 lb x 50 cm nylon suture
- 10 mm crimp

#### 2 x CCL Suture System (Swaged Needle)

- 45 mm fabella needle
- 80 lb x 80 cm nylon suture
- 12 mm crimp

#### 2 x CCL Suture System (Swaged Needle)

- 62 mm fabella needle
- 100 lb x 80 cm nylon suture
- 12 mm crimp



**PET-VET BIOMED Ltd.**  
32 Shaham St.  
Petach Tikva, Israel  
Tel: 972.3.9277290  
Fax: 972.3.9277278  
info@petvetBiomed.com  
www.petvetBiomed.com

**Show Room**  
8 Hamelacha St.  
Entrance A  
2nd floor, Suite 206  
South Industrial Area  
Netanya

**אולם תצוגה**  
רח' המלאכה 8  
כניסה א', קומה ב'  
חדר 206  
אזור תעשייה דרומי  
נתניה

**פוטו ביומד בע"מ**  
רח' שחם 32, קרית מטלון  
ת.ד. 7803, פתח תקוה 4951727  
טל.: 03.9277290  
פקס: 03.9277278  
info@petvetBiomed.com  
www.petvetBiomed.com