



CUTANPLAST

PET-VET BIOMED

Features

- ✓ Completely absorbed when left in situ
- ✓ Non-toxic
- ✓ Does not cause allergic reactions
- ✓ Has a neutral pH
- ✓ Absorbs blood corresponding to approximately 50 times its own weight
- ✓ May be applied dry
- ✓ May be soaked with sterile physiological solution or with antibiotics
- ✓ Is reabsorbed within a number of weeks
- ✓ Can easily be reduced in size to any required format
- ✓ Reduces the time of the surgical intervention
- ✓ Stops the hemorrhage
- ✓ Easy to use
- ✓ Practical
- ✓ Safe
- ✓ Reliable



PET-VET BIOMED Ltd.
32 Shaham St.
Petach Tikva, Israel
Tel: 972.3.9277290
Fax: 972.3.9277278
info@petvetBiomed.com
www.petvetBiomed.com

Show Room
8 Hamelacha St.
Entrance A
2nd floor, Suite 206
South Industrial Area
Netanya

אולם תצוגה
רח' המלאכה 8
כניסה א', קומה ב'
חדר 206
אזור תעשייה דרומי
נתניה

פט-וט ביומד בע"מ
רח' שחם 32, קרית מטלון
ת.ד. 7803, פתח תקוה 4951727
טל.: 03.9277290
פקס: 03.9277278
info@petvetBiomed.com
www.petvetBiomed.com