

For your safety during anaesthesia

OFFER OF
THE MONTH
July



21348801

EickSorber Anaesthetic Gas Absorber

Activated carbon filter to absorb organic anaesthetic gases

- ▶ Connects to pressure relief valve. Thus preventing unwanted gases from entering room atmosphere.
- ▶ Replace filters after exceeding 1,000 g.
- ▶ Includes 2 m exhaust hose

213488

EickSorber Single Unit

21348801



215205

Soda Lime Anaesthesia

Use with rebreathing circuits to prevent rebreathing of CO₂

- ▶ Low dust particles
- ▶ Colour indicator: white to violet
- ▶ 5 Litre

215205

**ORDER
ONLINE NOW**
www.eickemeyer.com