



QUICKCLEAN™ ULTRASONIC CLEANER

POWERFUL CLEANING – CONSISTENT RESULTS

Providing powerful, effective cleaning and consistent results, the QuickClean™ Ultrasonic Cleaner is a piece of equipment your surgery can't be without. The QuickClean's cutting-edge technology can help you create a safer, more efficient work environment by decreasing worker exposure to contaminants and potential sharps injuries while reducing the time and effort needed for cleaning.

LEAP

POWERFUL TECHNOLOGY

The QuickClean™ unit cleans thoroughly and consistently throughout the entire bath using Frequency-Leap technology. This advanced software helps ensure all of your instruments are fully cleaned the first time, every time, no matter where they are placed in the basket.



SIMPLE TO OPERATE

The QuickClean™ is easy to use right out of the box, so you can get your staff up and running with minimal training time. Features like a digital LCD control panel that displays prompts, options for preset or customized settings, intelligent software that remembers the last cycle settings, and a fluid sensor that ensures the bath always has sufficient water levels help simplify your processes and increase productivity.



RANGE OF SIZES AND DESIGN OPTIONS

The QuickClean™ comes in three table top sizes (4.5L, 12.5L, and 25L), so you can choose the option that best fits your configuration and workflow. And the stainless steel finish and attractive design coordinate well with any décor.

CUSTOMIZED
CYCLE TIMES FROM

1-99
MINUTES

**FULL-SIZE
SAFETY BASKET
INCLUDED**

**LIMITED
3-YEAR
WARRANTY***



Specifications, accessories and service part information on reverse.

* Full warranty terms & conditions available on request

SPECIFICATIONS:

QUICKCLEAN™ MODEL NUMBERS:	TABLE TOP UNIT: QC1-02	TABLE TOP UNIT: QC3-02	TABLE TOP UNIT: QC6-02
Capacity:	4.5 L	12.5 L	25L
Overall Length:	34.5 cm	37.5 cm	55 cm
Overall Width:	26 cm	42 cm	42 cm
Overall Height:	26.5 cm	26.5 cm	36.5 cm
Tank Size:	30 x 15 x 15 cm	33 x 30 x 15 cm	50.5 x 30 x 20 cm
Table Top Unit Weight:	5.76 kg (empty)	8.42 kg (empty)	13.44 kg (empty)
Table Top Shipping Weight:	7.5 kg	10.3 kg	15.6 kg
Operating Frequency:	32-38 kHz	32-38 kHz	32-38 kHz
Temperature up to:	60°C	60°C	60°C
Stainless Steel Safety Basket:	Included	Included	Included
Electrical Requirements:	Input: 230VAC +/- 10%, 50 Hz, 260W, 1.1A Output: 100W	Input: 230VAC +/- 10%, 50 Hz, 510W, 2.2A Output: 200W	Input: 230VAC +/- 10%, 50 Hz, 1060W, 4.6A Output: 400W
Standards:	BS EN 61010-1:2010 BS EN 60601-1-2:2007 BS EN 61326-1:2013 RoHS Compliant	BS EN 61010-1:2010 BS EN 60601-1-2:2007 BS EN 61326-1:2013 RoHS Compliant	BS EN 61010-1:2010 BS EN 60601-1-2:2007 BS EN 61326-1:2013 RoHS Compliant
Warranty:	3 Years	3 Years	3 Years

ACCESSORIES & SERVICE PARTS:



EXTRA SAFETY BASKETS

COM060100 – QC1 Safety Basket

COM060101 – QC3 Safety Basket

COM060102 – QC6 Safety Basket

Tel: +39 0522 875 166
newmed.sales@midmark.com
www.midmarkvetemea.com

Midmark EMEA Ltd, Beech House, First Floor, East Wing
Ancells Business Park, Fleet, Hampshire, GU51 2UN
Tel: +44(0) 1252 360 940 | hello@midmark.com

Manufactured by: Newmed Srl - a Midmark Company, Via Lenin 79/A,
42020 Quattro Castella, Reggio Emilia, Italy.
Tel: +39 0522 875 166 | Fax: +39 0522 243 096

UKM-121-00A1 0717

