

מוצר אורתופדי חדש של Eickemeyer

מחיר השקה – 2150 ₪ (במקום 2550 ₪)




PinCube

veterinary technology for life
Eickemeyer®

Precise Pin Positioning

NEW



PinCube

As part of an osteosynthesis procedure, bone and bone fragments are fixed for mutual stabilisation. This is usually done with plates, screws, tension bands and/or rods. For some applications, such as tension band wiring, it is recommended that k-wires are positioned parallel to each other. This is free hand hardly to be managed.

For this reason, EICKEMEYER® has developed the PinCube!

Advantages

- ▶ K-wires can be positioned exactly parallel
- ▶ The unique spiral-shaped hole pattern allows pins to be set both different distances and different diameters
- ▶ Usage of pins, drills and thread cutters in diameters of 0.6 – 2.0 mm
- ▶ Due to its extremely compact dimensions of maximum 35 mm, the PinCube is a very practical tool for small animal surgery

191348

הזמנות:

Natasha@petvetbiomed.com

shifra@petvetbiomed.com

T: 03-9277290 ;

F: 03-9277278



PET-VETBIOMED

veterinary technology for life
Eickemeyer®