

A warm nest

OFFER OF
THE MONTH
November



ICU for small patients

The Vetario T30 unit is an advanced, compact Intensive Care Unit capable of delivering small animal patient warming and recovery. It features full digital temperature control and early warning alarms. The Vetario T30 is ideal for young, sick, injured or weak animals in emergency situations or animal rescue. Its small profile makes it usable in a variety of veterinary settings, from pre-examination warming in consulting rooms, to use in patient homes or mobile ambulances.

- ▶ Digital temperature control, programmable to 0.1 °C
- ▶ Digital temperature reading – choose °C or °F
- ▶ High and low, internal and external temperature alarms
- ▶ Clear, high visibility locking-door
- ▶ Adjustable ventilation
- ▶ Made of anti-microbial plastics
- ▶ Easy to dismantle for thorough cleaning
- ▶ Built-in carry handle
- ▶ 12 V power supply for optional in-car adaptor
- ▶ Stackable with no assembly required
- ▶ Floor area: 240 x 240 mm
- ▶ Total volume: 9.5 litres
- ▶ Dimensions (in mm): H 297 x W 270 x D 275

219036

**ORDER
ONLINE NOW**
www.eickemeyer.com